
MAJOR
ACRYLIC
RESINS

major.base²⁰

**DENTURE BASE
POLY(METHYLMETHACRYLATE)-BASED
POLYMER FOR DENTAL PROSTHESIS.
HEAT-PROCESSED POLYMER.
POWDER AND LIQUID.**

major.base²⁰ has been developed in order to reach the highest odontotechnical standards and to exceed ISO 1567:2000 requirements (NIOM Certification of Conformity).

major.base²⁰ is available in the plain or veined version, in three different natural colours, selected to satisfy every clinical need.

major.base²⁰ performs high dimensional and colour stability and a great working ease. The fast-setting protocol drastically reduces waiting times during the polymerization procedure. In order to prevent allergic reactions to the product, **major.base²⁰** has a low tenor of residual monomer and it is well tolerated by patients.

INTENDED USE:

- Total prosthesis
- Partial prosthesis
- Clasp dentures

TECHNICAL STANDARD:


- ISO 1567:2000
- TYPE 1 (Heat-polymerized polymer)
- CLASS 1 (Powder and liquid)

POLYMERIZED MATERIAL TECHNICAL DATA:

- Flexural strength: 78 MPa
- Flexural modulus: 2390 MPa
- Water absorption: 22,0 µg / mm³
- Water solubility: 1,5 µg / mm³
- Residual monomer: 1,8 %

MAIN FEATURES:

- Available in a wide range of colours
- Fast polymerization protocol
- High dimensional stability
- High colour stability
- Low tenor of residual monomer

| AVAILABLE COLOURS | | |
|-------------------|------------------------|--|
| t | Transparent | |
| m | Pink / Red (plain) | |
| mv | Pink / Red (veined) | |
| i | Pink / Orange (plain) | |
| iv | Pink / Orange (veined) | |
| k | Dark Pink (plain) | |
| kv | Dark Pink (veined) | |

| SALE FORMAT: | | |
|--------------|---------------|--|
| code C3000 | Pack liquid: | 10 x 500 ml |
| code C3010 | Pack powder: | 10 x 500 g |
| code C3020 | Intro set: | Powder 100 g Liquid 100 ml |
| code C3021 | Standard set: | Powder 2 x 500 g Liquid 1 x 500 ml Plastic box |

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major.repair

**POLY(METHYLMETHACRYLATE)-BASED
POLYMER FOR DENTAL PROSTHESIS.
AUTOPOLYMERIZABLE POLYMER.
POWDER AND LIQUID.**

major.repair has been developed in order to reach the highest odontotechnical standards and to exceed ISO 1567:2000 requirements (NIOM Certification of Conformity).

major.repair performs high dimensional and colour stability, a great working ease and a perfect chemical bonding to denture base polymers.

major.repair is available in the plain or veined version, in three different natural colors, selected to match major.base₂₀ shades.

major.repair can be polymerized as a pour-type resin or under pressure, giving excellent results in both techniques.

In order to prevent allergic reactions to the product, **major.repair** has a low tenor of residual monomer and it is well tolerated by patients.

INTENDED USE:

- To repair and rebase prosthesis
- Temporary prosthesis

TECHNICAL STANDARD:

- ISO 1567:2000
- TYPE 2 (Autopolymerizable polymer)
- CLASS 1 (Powder and liquid)

POLYMERIZED MATERIAL TECHNICAL DATA:

- Flexural strength: 66,4 MPa
- Flexural modulus: 2217 MPa
- Water absorption: 21,2 µg / mm³
- Water solubility: 1,4 µg / mm³
- Residual monomer: 4,0 %

MAIN FEATURES:

- Colours matching Major.Base₂₀ shades
- High dimensional stability
- High colour stability
- Low tenor of residual monomer

| AVAILABLE COLOURS | | |
|-------------------|------------------------|--|
| m | Pink / Red (plain) | |
| mv | Pink / Red (veined) | |
| i | Pink / Orange (plain) | |
| iv | Pink / Orange (veined) | |
| k | Dark Pink (plain) | |
| kv | Dark Pink (veined) | |

| SALE FORMAT: | | |
|--------------|---------------|--|
| code C3100 | Pack liquid: | 10 x 500 ml |
| code C3110 | Pack powder: | 10 x 500 g |
| code C3120 | Intro set: | Powder 100 g Liquid 100 ml |
| code C3121 | Standard set: | Powder 2 x 500 g Liquid 1 x 500 ml Plastic box |

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major.ortho

POLY(METHYLMETHACRYLATE)-BASED POLYMER FOR THE CONSTRUCTION OF REMOVABLE ORTHODONTIC APPLIANCES. AUTOPOLYMERIZABLE POLYMER. POWDER AND LIQUID.

major.ortho performs high dimensional stability and a great working ease.

major.ortho has been developed to be employed in three different techniques, always assuring the best results:

- Major.Ortho can be polymerized under pressure for the construction of bites and removable orthodontic appliances.
- Major.Ortho is suitable for the "pepper and salt" technique for the construction of orthodontic plates or rehabilitative bites.
- Major.Ortho can be applied using a brush to realize small parts or to repair orthodontic appliances.

INTENDED USE:


- Removable orthodontic appliances
- Prosthetic "bites"
- Transparent plates for teeth mounting

POLYMERIZED MATERIAL TECHNICAL DATA:


- Flexural strength: 50 MPa
- Flexural modulus: 1770 MPa
- Water absorption: 17,0 µg / mm³
- Water solubility: 3,6 µg / mm³

MAIN FEATURES:

- High dimensional stability
- Suitable for different working techniques
- NIOM certified

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major.tray

**POLY(METHYLMETHACRYLATE)-BASED
POLYMER FOR THE CONSTRUCTION OF
INDIVIDUAL IMPRESSION TRAYS.**

**AUTOPOLYMERIZABLE POLYMER.
POWDER AND LIQUID.**

major.tray performs high dimensional stability and a great working ease.

It can be immediately manipulated after mixing and it is bench self-polymerizable (in ten minutes after mixing).

major.tray is available in three different colours (yellow, green and blue).

INTENDED USE:

- Individual impression trays
- Various supports

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- Breaking strength: 21 N
- Water absorption: 85,0 µg / mm³
- Water solubility: 7,5 µg / mm³

MAIN FEATURES:

- Available in a wide range of colours
- Bench self-polymerization
- High dimensional stability
- Working ease

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MAJOR PRODOTTI DENTARI S.p.A. - via L. Einaudi, 23 - 10024 Moncalieri (To)